



Mains pressure indoor stainless steel coiled cylinders

This diagram is provided to assist identifying parts if replacement is necessary. Contact your place of purchase or the manufacturer to enquire about parts, availability and/or service.


Element replacement

The coiled cylinders had an element change in July 2021. To ensure the correct element is ordered check the date of manufacture on the cylinder label. The cylinder may also have a spares label attached, refer below image. If it does, order the code shown on the label.



ELEMENT SPARE PARTS LIST	
PART NO	DESCRIPTION
30060	ENAMEL 2KW
30061	ENAMEL 3KW
30062	STAINLESS STEEL 2KW
30063	STAINLESS STEEL 3KW
(* includes element seal & backnut)	
Torque tighten backnut to 12Nm	

30062: Kawai element 2 kW
(has a temp. sensor pocket)



For cylinders made Jul 2021 onwards, includes seal and backnut.

30063: Kawai element 3 kW
(has a temp. sensor pocket)



For cylinders made Jul 2021 onwards, includes seal and backnut.

Note: If you damage the Kawai element backnut and/or seal and need a replacement, these can be ordered separately.



20936: Element backnut Kawai



20937: Element seal Kawai

30002: 2 kW curved element



For cylinders made **BEFORE** July 2021, includes seal and backnut.

30004: 3 kW curved element



For cylinders made **BEFORE** July 2021, includes seal and backnut.

30006: 3 kW curved titanium element



Alternate element option for hard and aggressive water areas.

20356 Backnut element, to replace older or damaged nuts for 30002 / 30004



30022 Element gasket (rubber seal), to replace leaking or damaged seals for 30002 / 30004. Two seals provided in a kit.



30015: Thermostat Robert Shaw orange dial (ART)




30018: TPR mains pressure 850 kPa 15 mm



For MP single **non-thermosiphon** cylinders.

BE11004781
TPR MP 850 kPa 20 mm



For MP single and twin coil **thermosiphon** cylinders

30019: Element cover box



To replace old or damaged covers.

20355: Blanking plate incl. gasket



If the top or bottom element is not needed the blanking plate can be used to 'blank' the unused connection.

30020: Element key (to remove element backnuts)



As of June 2023, the element key is no longer included with the cylinder; it is now part of the element spare kit. DO NOT leave the element key inside the element enclosure due to potential electrical hazards.